

## GEWA-FIT

- GEWA-FIT claw and hoof treatment set is ideal to calm injured or diseased claws
- healing procedure of hoof diseases such as sole ulcers, defects of the corium, abscesses on soles, etc., is supported in a positive way
- the special glue is skin-friendly, toxically harmless and easy to apply. Ideal treatment temperature between 15°C and 20°C
- the glue hardens between 5 - 7 minutes
- the animals get accustomed to the wooden block easily, even lame cows are able to walk much better after short time
- GEWA-FIT treatment sets are available in packs of 5, 10 or 25 pcs.
- choosing GEWA-FIT you receive a high-quality product, which has been tested to German DLG standards
- optimal price-performance ratio
- GEWA-FIT supports the vet's therapy and shortens the often quite long treatment times significantly

## GEWA-FIT treatment set, complete with 25 wooden blocs

- sufficient for 25 - 35 treatments → contents: 1 kg of powder, 1 l of fluid
- 25 wooden blocs, 25 spatula, 1 mixing cup, 1 dosing spoon



## GEWA-FIT treatment set, complete with 10 wooden blocs

- for 10 to 14 treatments → contents: 500 gram of powder, 500 ml of fluid →
- 10 wooden blocs, 10 spatula, 1 mixing cup, 1 dosing spoon



## GEWA-FIT treatment set, complete with 5 wooden blocs

- for 5 - 7 treatments → contents: 250 gram of powder, 500 ml of fluid
- 5 wooden blocs, 5 spatula, 1 mixing cup and 1 dosing spoon



Type	Ref.		
GEWA-FIT treatment set, complete with 25 wooden blocs	300.0301		
GEWA-FIT treatment set, complete with 10 wooden blocs	300.0306		
GEWA-FIT treatment set, complete with 5 wooden blocs	300.0311		

